Comfort Cool® Pronation-Supination Splint

Instructions for Care

Hand wash in cool water with mild soap and air dry. Fasten strap during washing.

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
Left	Right	Size	MP Circumference
NC34800	NC34801	Pediatric	5½" to 6¼" (14 to 16 cm)
NC34802	NC34803	Small	6¼" to 7" (16 to 18 cm)
NC34804	NC34805	Medium	7" to 8" (18 to 20 cm)
NC34806	NC34807	Large	8" to 9" (20 to 23 cm)

MD

60°F (16°C) — 85°F (29°C)

North Coast Medical, Inc.
135 E Main Ave, Suite 110
Morgan Hill, CA 95037 – U.S.A.

EC REP

ICON (LR) Limited
South County Business Park
Leopardstown, Dublin 18
D18 XSR3, Ireland

UK REP

ICON (LR) Limited
500 South Oak Way
Green Park, Reading
RG2 6AD, United Kingdom

CH REP

Allenspach Medical AG
Dornacherstrasse 5
4710 Balsthal, Switzerland

Comfort Cool® Pronation-Supination Splint

Lightweight, comfortable support.

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Attelle de pronation-supination

Pronations-/Supinations-Schiene

Pronations-/supinationssskena

Helps position the upper extremity for improved function, supports thumb prehension, and reduces abnormal tone.

The strap provides dynamic support for supination or pronation of the forearm.

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Comfort Cool® Pronation-Supination Splint

Intended Purpose

Positions the upper extremity for better function, facilitates prehensive use of the thumb, and reduces abnormal tone.

Indications

Suitable for neurological conditions such as CP, CVA, and TBI, brachial plexus injuries, and orthopedic conditions.

Effective for opening a shortened web space due to soft tissue injuries, thumb tendinitis, nerve palsy, or contracture from disuse.

Contraindications

Not for use with severe spasticity, significant fluctuations in tone, or edema.

Do not use if allergic to neoprene or synthetic fibers.

Instructions for Use

Hand Piece:

1. Slide the thumb and index finger into the hand piece with the web space snug up against the inside seam.
2. Wrap the wrist strap around to the back of the hand. To adjust thumb abduction, tighten or loosen the radial side of the hand piece over the wrist strap. (I)
3. Trim the neoprene for a custom fit, as needed. The seams will not unravel.

Achieving Supination:

4. Place the extended hook end of the forearm strap diagonally on the palm side of the hand piece. (II)
5. Pull the strap through the thumb webspace toward the dorsal hand. Attach the strap's second hook piece onto the back of the hand piece. (III) If the web strap is too wide for thumb comfort, trim the edges of the neoprene strap to make it narrower.
6. Wrap the strap around the wrist on the ulnar side, spiraling up to the elbow. Wrap the end of the strap around the elbow area at least once to form a cuff, securing with the separate hook piece. This cuff may be placed just above or just below the elbow, as needed. (IV)
7. To create a stronger pull, hold the arm in the desired supinated position. Rewrap the spiral strap, starting at the hand, pulling it taut as you progress up to the elbow. Trim if needed and refasten the cuff.
8. If alternating between pronation and supination, trim the strap based on supination length to avoid over-trimming.

Comfort Cool® Pronation-Supination Splint

Instructions for Use

Achieving Pronation:

1. Place the extended hook end of the forearm strap diagonally onto the back of the hand piece. (V)
2. Pull the strap through the thumb web space, toward the palm. Fasten the strap's second hook piece onto the hand piece. If the web strap fits too wide / the thumb is uncomfortable, the edges of the neoprene strap can be trimmed. (VI)
3. Continue with Supination steps #6-8 (at left) but wrap in the opposite direction to achieve pronation. (VII)

If no special thumb position is desired:

If the thumb does not need to be placed in more abduction or opposition, the end hook section of the spiral forearm strap can be cut off at the second hook piece on the strap. Attach this second hook piece diagonally on the hand piece.

For **supination**, place it diagonally over the dorsal wrist area, pulling the strap toward the ulnar side. (VIII)

For **pronation**, place it diagonally over the palmar area, pulling the strap toward the ulnar side. Wrap it spirally up the arm. The straps do not pass through the web space with this design.