



Intended Purposes Helps prevent and/or limit ulnar drift of the fingers. Provides gentle support along the palmar MP joints while the individual extension pieces dynamically pull the proximal phalanx's radially into better alignment. Fabricated with loop material on the outside for easy fastening. The thick neoprene design is warm and holds fingers more firmly.

Indications Rheumatoid arthritis, degenerative joint disease.

Contraindications Skin irritation, allergy to nylon, neoprene, polyester, or spandex, loss of circulation.

Instructions for Use / Proper Fit



Place the wide middle section of the support across the palm with the end (with sizing label) passing through the thumb web space. Fasten both ends together on the back of the hand (I). The support should be just proximal to the MP joint/knuckles (Figure 2). Position the first palmer extension piece (II) so it lines up between the index and middle finger (Figure 1). Bring the extension piece through the index's web space, then around the index finger and pull it diagonally toward the radial edge (thumb) of the hand. Fasten it in place on the wide support (III) with the hook tab. Continue in the same fashion, wrapping the next extension piece around the next finger and attaching it to the previous extension, overlapping each in a radial direction. The last extension piece wraps across the top of the little finger to attach to the previous extension (IV). Splint may be trimmed as needed. Avoid cutting into stitching.

Warning: If extension pieces feel too tight, affecting circulation, loosen as needed. If fingers change color, there is an increase in swelling, skin becomes irritated, or if pain, stiffness, or numbness increases, discontinue use and seek advice from a healthcare professional.

Instructions for Care Hand wash in cool water and mild soap. Air dry. Do not use dryer. Machine washing and drying can damage the support.