

### Intended Purposes

Non-sterile device used to assist in pulling on sock or stockings, with use of a single or double cord.

### Indications

Persons with difficulty bending over, dizziness, or vertigo.

### Contraindications

Limited hand strength. Allergy to polypropylene.

### Instructions for Use



1 Gather entire sock over end of sock aid.



2 Hold strap and drop sock aid onto floor.



3 Slide foot into sock aid opening.



4 Pull upward on strap to bring sock over heel and up leg.

### Warning:

Keep out of reach of children. Plastic bag and small parts may pose a choking hazard.

### Instructions for Care

Hand wash in cool water and mild soap. Air dry.

*Note: Cord handle can be cut and securely knotted to adjust length.*

