



# North Coast

MEDICAL & REHABILITATION PRODUCTS

EXPANDING HORIZONS IN  
PHYSICAL THERAPY

APTA CSM  
HOUSTON

FEBRUARY 13-15, 2025





YOUR **ONE** BRAND FOR  
DRY NEEDLING!

**BLUE DIAMOND™**  
DRY NEEDLES



**8 Needle Lengths / 9 Needle Diameters / 21 Needle Combinations**

- Features an oversized, conductive aluminum etched handle, ensuring superior control when utilizing enhanced target-specific techniques.
- Japanese high-grade surgical steel needle body. Designed to allow maximum flexion and deeper penetration.
- Medical-grade silicone lubricated needle body engineered to glide effortlessly and reduce friction – remarkably gentle for improved patient comfort.
- Each box contains 20 sleeves of 5 needles (100 pieces per box).
- FDA cleared.

*Easy to peel  
medical dialysis  
paper backing*



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Oversized, conductive, etched aluminum handle




Double-clamped tube for easy one-handed release



**BLUE DIAMOND™**  
DRY NEEDLES

## NEEDLE DIAMETER COLOR CHART



### 0.16 to 0.22mm Diameter

	Diameter	Length
NCBD-1613	0.16mm	13mm
NCBD-1813	0.18mm	13mm
NCBD-2013	0.20mm	13mm
NCBD-2213	0.22mm	13mm

### 0.25mm Diameter

	Diameter	Length
NCBD-2513	0.25mm	13mm
NCBD-2530	0.25mm	30mm
NCBD-2540	0.25mm	40mm
NCBD-2550	0.25mm	50mm
NCBD-2560	0.25mm	60mm
NCBD-2575	0.25mm	75mm
NCBD-25100	0.25mm	100mm
NCBD-25125	0.25mm	125mm


### 0.30 to 0.35mm Diameter

	Diameter	Length
NCBD-3030	0.30mm	30mm
NCBD-3040	0.30mm	40mm
NCBD-3050	0.30mm	50mm
NCBD-3060	0.30mm	60mm
NCBD-3075	0.30mm	75mm
NCBD-30100	0.30mm	100mm
NCBD-3575	0.35mm	75mm

### 0.40 to 0.45mm Diameter

	Diameter	Length
NCBD-40100	0.40mm	100mm
NCBD-45125	0.45mm	125mm

Reserved for licensed medical professionals and approved resellers only.





# Rainbow™

LATEX FREE  
EXERCISE TUBING

**Consistent Resistance  
Poundage Force**  
*while stretching!*



LEVEL 1

LEVEL 2

LEVEL 3

LEVEL 4

LEVEL 5



Color-Coded Resistance Levels with Functional Colors™

- No Residue
- No Powder
- No Odor



## Rainbow™ Latex-Free\* Exercise Tubing

### Professional

#### 100 ft Box

One 100 ft (30m) roll

NC95421-100	Level 1 - Light	Orange
NC95422-100	Level 2 - Medium	Pink
NC95423-100	Level 3 - Heavy	Green
NC95424-100	Level 4 - Extra Heavy	Blue
NC95425-100	Level 5 - Super Heavy	Purple

### Retail

#### Multipacks

Three 5 ft (1.5 m) lengths of tubing per pack

NC95430	Light Resistance Pack	Orange, Pink, Green	Levels 1, 2 and 3
NC95431	Heavy Resistance Pack	Green, Blue, Purple	Levels 3, 4 and 5



\*Not made with natural rubber latex.

To order, log in at  
**www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™



# Rainbow™

LATEX FREE  
EXERCISE BAND

A premium  
**Latex Free\* Exercise Band**  
without a premium price!

LEVEL 1

LEVEL 2

LEVEL 3

LEVEL 4

LEVEL 5



Color-Coded Resistance Levels with Functional Colors™

- No-Slip Grip
- No Powder
- No Odor



## Rainbow™ Latex-Free\* Exercise Band

### Professional

#### 50 yd Box

One 50 yd (45 m) roll

NC95401-50	Level 1 - Light	Orange
NC95402-50	Level 2 - Medium	Pink
NC95403-50	Level 3 - Heavy	Green
NC95404-50	Level 4 - Extra Heavy	Blue
NC95405-50	Level 5 - Super Heavy	Purple

### Retail

#### Multipacks - Three 5 ft (1.5m) lengths of band per pack

NC95410	Light Resistance Pack	Orange, Pink, Green	Levels 1, 2 and 3
NC95411	Heavy Resistance Pack	Green, Blue, Purple	Levels 3, 4 and 5

#### Complete Multipack - Five 5 ft (1.5m) lengths of band per pack

NC95400	Orange, Pink, Green, Blue, Purple	Levels 1, 2, 3, 4 and 5
---------	-----------------------------------	-------------------------



\*Not made with natural rubber latex.



To order, log in at  
**www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™



# Rainbow™ Exercise Band Mini Loops

Short, continuous band loops are ideal for upper and lower body exercises.




- Rainbow™ Exercise Band Mini Loops are seam-free and flat.
- Mini Loops measure 2" (5 cm) x 11" (28 cm).
- Available in five color-coded resistances.
- Contains latex.
- Sold individually.



## Rainbow™ Exercise Band Mini Loops

	Resistance	Color
NC91241-01	Level 1	Orange
NC91242-01	Level 2	Pink
NC91243-01	Level 3	Green
NC91244-01	Level 4	Blue
NC91245-01	Level 5	Purple



To order, log in at [www.ncmedical.com](http://www.ncmedical.com)  
or call **800.821.9319**

Available in Virtual Showroom™



# Exercise Balls

Economical burst-resistant balls provide versatile exercise workouts.

- Tested to support up to 550 lbs. (250 kg).
- Helps to improve core stabilization by training the spinal column muscles as the body works to stay in balance.
- Balls are available individually in retail box with handle.
- Not intended for heavy lifting.
- Not made with natural rubber latex.

## Exercise Balls

	Color	Diameter
NC50100	Blue	45 cm (17¾")
NC50101	Red	55 cm (21½")
NC50102	Teal	65 cm (25½")
NC50103	Dark Blue	75 cm (29½")
NC50104	Silver	85 cm (33½")



*Ideal for clinic or home use.*



# Mini Ball

Strengthen key muscle groups, core health, endurance, posture and balance.

- High quality, 9" (23 cm) diameter therapy ball intended for exercising and strengthening the abdomen, back, hips, knees, shoulders, and neck.
- Anti-burst vinyl material deflates slowly if punctured.
- Designed to remain elastic and squishy when fully inflated.
- Supports up to 330 lbs. (150 kg).
- Not made with natural rubber latex.



*Includes exercise poster.*

NC50109 Mini Ball

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™






# Economy Shoulder Pulley

Highest quality in the market!

- Increase shoulder range of motion in the clinic or at home.
- Large, easy-grip handles.
- Non-slip web door strap.
- Simple to adjust cord length.
- Smooth action pulley mechanism.



Economy  
Shoulder Pulley  
NC52080



# Universal Door Bracket

Easy on, easy off  
over-the-door bracket.

- Bracket can remain in place on the top of the door, ready for the next exercise.
- Made of smooth, coated metal to help prevent door scratches.
- Bracket width is 1<sup>3</sup>/<sub>4</sub>" (4 cm), which fits over commercial and most standard doors.
- Includes two EVA foam pads which help secure the bracket around doors with narrow widths.

NC52081 Universal Door Bracket



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™

# Shoulder Pulley

Everything needed to perform range of motion exercises at home.

- Ideal for home use, to independently improve range of shoulder movement and coordination following shoulder injury or surgery.
- Can be used for exercises while seated, standing, or lying down.
- Sets up in minutes by positioning the door strap over the top along the hinged side of any standard door.
- Pre-assembled and ready to use.
- The door strap fits snugly between the door and frame and will not scratch surfaces.
- For optimal attachment to door, use the Norco<sup>®</sup> Universal Door Bracket (sold separately).
- Not made with natural rubber latex.



NC52069 Shoulder Pulley



### **Retail ready box includes:**

- Exercise booklet with step-by-step instructions and photos
- Door strap
- 6½ ft. (2.4 m) cord
- Pulley mechanism and comfortable, ergonomically designed handles

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™



# Stretch Strap

Safe and effective deep, gradual stretching.



**ISOTONICS**  
**ISOMETRICS**  
**PROLONGED STRETCH**  
**INCREASED R.O.M.**  
**FLEXIBILITY**

- Designed with ten multi-positioned grip loops for use in deep, gradual stretching of the major muscle groups.
- Ideal for individuals seeking unassisted or assisted stretching exercises to relieve sore or achy muscles.
- Stretch more safely and effectively than is possible unaided.
- Grip loops can be used to set stretch improvement goals.
- Made of non-stretch webbing.
- Machine washable.
- Not made with natural rubber latex.

*Retail ready packaging includes an instruction booklet for quick reference.*



## Stretch Strap

	Model	Color	Length
NC84527	Standard	Purple	6' 4" (1.9 m) for people under 6' (1.8 m)
NC84528	Long	Blue	7' 5" (2.3 m) for people over 6' (1.8 m)

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™



Common Trigger Points



# BODYtool™

TRIGGER POINT SELF MASSAGER

## Repair Yourself

- Designed to target trigger points and muscle knots in the body to assist with temporary release of muscular tension.
- **Textured grip handles** allow better control for effectively targeting trigger points.
- Features **seven large knobs**, strategically placed to control pressure to the desired areas for trigger point release and breaking up of myofascial tissue.



Easy-to-follow instruction guide and poster included!

**NC69700 Body Tool™**

**Body Tool™  
Retail Stand  
NC69700-RS**



Retail display rack holds 12 Body Tool™ massagers.

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™




# EVA Foam Roller

**Professional quality foam roller designed for heavy and repeated use in multi-user environments.**

- Great for proprioceptive balance and stability exercise or to roll out trigger points, knots and areas of increased muscle density.
- Use for strengthening, flexibility and stretching exercises.
- Provides comfortable support — gentle on tight muscles while lying supine, yet functional enough for light massage rolling.
- Closed-cell, textured, high-density EVA foam technology resists breaking down and allows for easy cleaning between clients.
- Built with high-grade materials for long-lasting performance.
- Roll measures 36" (91 cm) long with a 6" (15 cm) diameter.

**EVA Foam Roller**  
NC84800-01



Not made with  
natural rubber latex.

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™



# Lumbar Roll


**Provides comfortable support to help improve posture and relieve lower back pain.**

- Adjustable, non-elastic strap allows users to secure the Lumbar Roll to standard chair backs and car seats at the desired height.
- Innovative strapping design allows the Lumbar Roll to maintain position in the hollow of the back, or be repositioned for maximum comfort.
- Made of polyurethane foam with a removable polyester cover that is easy to clean.
- The removable cover includes a zipper guard for easy opening and closing of the zipper.
- Available in **Standard** and **Firm** densities.
- Rolls measure 11" (28 cm) long with a 5" (13 cm) diameter.
- Not made with natural rubber latex.



## Lumbar Roll

	Foam Density
NC89095	Standard
NC89096	Firm



*Strap "tunnel" allows the strap to be tightened firmly to maintain correct positioning.*



**SCAN**  
to learn more!



*Zipper guard for snag-free zipper use.*

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™





# MUSCLE AIDTAPE™

**ARM YOUR MUSCLE**

**BEAT THE PAIN**

- Relieve muscle pain, soreness and cramping.
- Provide support for joints and muscles.
- Reduce swelling, inflammation and bruising.
- Increase blood circulation/flow.
- Water-resistant and breathable.
- Not made with natural rubber latex.



## Available in:

Single rolls: 2" x 16' (5 cm x 5 m)

3" x 16' (10 cm x 5 m)

4" x 16' (10 cm x 5 m)

Bulk rolls: 2" x 105' (5 cm x 32 m)

Display Cartons	Color	Quantity	Width
MT1642-MVBG	Beige	6 Rolls	2" (5 cm)
MT1642-JDBK	Black	6 Rolls	2" (5 cm)
MT1642-RPOR	Orange	6 Rolls	2" (5 cm)
MT1642-MVPK	Pink	6 Rolls	2" (5 cm)
MT1642-JDNB	Navy Blue	6 Rolls	2" (5 cm)
MT1642-JDBL	Blue	6 Rolls	2" (5 cm)
MT1642-RPUS	USA	6 Rolls	2" (5 cm)
MT164W3-BK	Black	4 Rolls	3" (8 cm)
MT164W4-BK	Black	3 Rolls	4" (10 cm)

Bulk Rolls	Color
MT1052-BG-BULK	Beige
MT1052-BK-BULK	Black
MT1052-BL-BULK	Blue



*Also Available in*  
**10" Long Precut Strips**

**Strengthen Performance.  
 Relieve Muscle Pain.  
 Facilitate Rapid Recovery.**



Precut strips: 2" x 10" (5 cm x 25.5 cm)

<b>Precuts (20)</b>	<b>Color</b>	<b>Bulk Precuts (128)</b>	<b>Color</b>
MTPCBG-20	Beige	MTPCBG-128	Beige
MTPCBK-20	Black	MTPCBK-128	Black
MTPCBL-20	Blue	MTPCUS-128	USA
MTPCUS-20	USA		



**Reid Priddy**  
 Olympic Gold Medalist

**MUSCLE  
 AIDTAPE™**

*Longer lasting, skin-friendly support.*




To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

# Inclinometer

## Brake button holds range of motion results with ease.

- For use in measuring joint angles, including range of motion (ROM) and movement.
- Facilitates assessment and monitoring during physical rehabilitation or medical treatment.
- Ideal for measuring flexion, extension, abduction, adduction, as well as rotation in the neck, shoulder, elbow, wrist, hip, ankle and spine.
- Adjustable dial measures range of motion traveled.
- Instruction manual includes 22 exercise samples.




**NC70309** Inclinometer



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**





# MAXIMUM PERFORMANCE PAIN RELIEF.

## CLINICAL

- Provides fast, long-lasting relief from sprains, strains, bruises, arthritis, muscle aches, joint and back pain.
- Specially formulated penetration enhancers to provide fast pain relief.
- Contains 10% menthol, dimethyl sulfone (MSM) and glucosamine.
- Sweat-resistant.
- Paraben-free, no dyes or fragrances.
- Non-greasy and non-staining.



### Stopain® Clinical - Singles Qty

NC80091-02	Roll-On	3 fl. oz. (88 ml)	1
NC80091-03	Spray	4 fl. oz. (118 ml)	1
NC80091-01	Gel Tube	4 fl. oz. (118 ml)	1

### Stopain® Clinical - Larger Size Gel Pump Qty

NC80091-05	32 fl. oz. (946 ml)	1
NC80091-06	1 gallon (3.78 liters)	1
NC80091-04	16 fl. oz. (473 ml)	1

### Stopain® Clinical - Starter Pack Qty

NC80091-08	Sample Packs w/ Display Box	(100) 3 gram packets
------------	-----------------------------	----------------------



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™



# Eco-Stim™ TENS + EMS

Combo unit provides safe, drug-free pain relief.

- TENS, EMS and Massage wave forms.
- 22 treatment modes. 6 modes allow Time, Frequency and Pulse Width customization.
- Large, easy-to-read backlit display.
- 40 levels of output intensity.
- Safety feature locks intensity level.
- Dual channel 120mA (at 500 Ohm load).
- Output frequency:  
2 - 100 Hz  
Pulse width:  
100- 300  $\mu$ s per pulse
- Includes:  
(4) 2" x 2" square reusable electrodes,  
(2) lead wires,  
(3) AAA batteries,  
(1) belt clip, and  
(1) instruction manual.

NC89481  
Eco-Stim™ TENS + EMS

*Use belt clip for  
on-the-go pain relief.*

Not made with  
natural rubber latex.

*Includes 48 page  
instruction manual.*



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™




Latex Free\*

# Eco-Stim™

## Neurostimulation Electrodes



- Premium quality Carbon/Hydrogel for maximum comfort.
- Reusable up to 20 times when stored as recommended.
- Compatible with all NMS and TENS devices.
- Personalized packaging for patient data usage.
- FDA cleared for OTC purchase.
- 2mm DIN pigtail connector.
- Available in a variety of shapes and sizes in both foam or cloth
- Not made with natural rubber latex.




**2" x 2"  
Square**

Foam: NC89292-1  
Cloth: NC89293-1




**2" x 4"  
Rectangle**

Foam: NC89292-3  
Cloth: NC89293-3



**2" diameter  
Round**

Foam: NC89292-2  
Cloth: NC89293-2



**3/4" diameter  
Round**

Cloth: NC89293-4

Cloth

Foam

*Personalized packaging  
for patient data usage.*

*4 electrodes  
per pack*



To order, log in at  
**www.ncmedical.com**  
or call **800.821.9319**  
Available in Virtual Showroom™

# Wireless "On-the-Go" Technology

Intelligent iontophoresis.



Technologically advanced, self-contained iontophoresis drug delivery system.



- Fully self-contained.
- Intelligent biofeedback microprocessor.
- Hydrogel dispersive technology—no saline required.
- Patented "shut-off" technology.
- SmartPower™ LED.
- Detects and measures skin impedance to optimize drug delivery.



**Active area:** 11.4 cm<sup>2</sup>  
**Dosage:** 40 or 80 mA\*min  
**Skin Adhesive:** Hypoallergenic  
**Maximum Buffer Capacity:** Ag/AgCl Stabilizer  
**Fill Volume:** 2.0cc  
**Smart Power:** 6V with intelligent microprocessor  
**Packaging:** 6 ActivaPatch® Kits, Alcohol Wipes, Patient Information Cards  
**Treatment time:** 2.5 Hours

ActivaPatch® IntelliDose™ 2.5  
NC89188

**80 mA\*min dose  
in only 2.5 hours.**

\* Study data on file

**no saline needed.**

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319



## ActivaPatch® IontoGo™ technology.

No external power source to connect, no leadwires,  
no dose controller needed, power shut down feature and  
**no saline required.**

**ACTIVApatch®**  
IontoGo™ 4.0

Extended wear patch with an average drug delivery dose of 80 mA\*min in approximately 4 hours.

**4**  
Hours

**ACTIVApatch®**  
IontoGo™ 12.0

Extended wear patch with an average drug delivery dose of 80 mA\*min in approximately 12 hours.

**12**  
Hours

**HOW SMART  
IS YOUR  
WIRELESS IONTO?**

**Active area:** 11.4 cm<sup>2</sup>  
**Dosage:** 80 mA\*min  
**Skin Adhesive:** Hypoallergenic  
**Fill Volume:** 2.0cc  
**Packaging:** 6 IontoGo™ 4.0,  
6 Alcohol Wipes  
**Treatment Time:** 4 Hours

**NO  
SALINE  
REQUIRED**

**ActivaPatch®  
IontoGo™ 4.0  
NC89189**

**FULLY  
SELF-  
CONTAINED**



**Active area:** 11.4 cm<sup>2</sup>  
**Dosage:** 80 mA\*min  
**Skin Adhesive:** Hypoallergenic  
**Fill Volume:** 2.0cc  
**Packaging:** 6 IontoGo™ 12.0,  
6 Alcohol Wipes  
**Treatment Time:** 12 Hours

**NO  
SALINE  
REQUIRED**

**ActivaPatch®  
IontoGo™ 12.0  
NC89191**

**FULLY  
SELF-  
CONTAINED**



\* Study data on file

# Smarter Iontophoresis™

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

# Wired Technology



## ADVANTAGES

- Conforms extremely well to uneven treatment sites.
- Maximum strength, precise buffering with new Carbon-Ag/AgCl technology.
- Large ground electrode surface area lowers sensation and increases comfort.
- Highly absorbant drug matrix.
- Large surface adhesion area.



## TRIVARION® IS AVAILABLE IN THE FOLLOWING SIZES:



SMALL 1.5cc  
NC89180B



MEDIUM 2.0cc  
NC89181B



LARGE 3.5cc  
NC89182B

TRIVARION®  
BUTTERFLY



MEDIUM 2.0cc  
NC89183B



ActivaDose® II Controller  
NC89186

## ADVANTAGES

- Two controls operate dose and current, with optional push-knob or turn capabilities.
- High-contrast display indicates dose, time and current.
- More power than dual-channel devices for better penetration; up to 80 volts.
- Automatic current ramp-up and ramp-down for maximum comfort. Dose range from 0 to 80 mA-min with 0.1 mA-min increments.
- Time calculations are performed automatically and recalculated if dose and current settings are changed, or with pause/restart.

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319



# Buffered Iontophoretic Delivery Kit



Each electrode has a large surface adhesion area.

Variety of sizes and shapes offer great flexibility for iontophoresis treatments.

- Made of non-sensitizing, hypoallergenic, medical grade materials.
- Delivers ions with a pH balance up to 80 mA-min.
- Maximum strength, precise buffering with unique carbon- Ag/AgCl technology.
- The highly absorbent drug matrix has no messy hydrogel.
- Large ground electrode surface area lowers sensation and increases comfort.
- Each kit contains 12 individual treatment sessions that include: one activity electrode; one return electrode; one alcohol prep pad.

Flexible. Thin. Contouring. Comfortable.



Competitor's product is not as thin or flexible.



Advanced, technology-based manufacturing for a thinner, more pliable electrode.



Easy-to-open packaging.

Not made with natural rubber latex.



**Small**  
Square 1.5 cc  
NC89250B



**Medium**  
Square 2.5 cc  
NC89251B



**Large**  
Square 4.0 cc  
NC89252B



**Butterfly 2.0 cc**  
NC89255B

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**



# 3-Section Hi-Lo Table Electric

Full-feature table for ease of use and client comfort.

- Dual full-length foot bars and hand controller for easy height adjustment.
- Premium matte black polyurethane upholstery offers comfort and durability. Easy to clean and sanitize.
- Superior comfort with 2.5" foam padding. (1.25" standard foam, plus 1.25" Slo-Foam™)
- Removable face plug with contoured opening.
- Convenient paper roll dispenser and cutter for quick setup between clients.
- Durable white powder-coated steel frame.
- One-step caster control bar engages all casters at once for easy table repositioning.

**3-Section Hi-Lo Table, Electric**  
NC94710



## Specifications:

- 30" (76cm) Wide; 76" (193cm) Length
- Height Range 19" to 37" (48cm to 94cm)
- Lifting Capacity: 550 lbs (250 kgs)
- Static Weight Capacity: 850 lbs (385 kgs)
- Head section angle ranges from -90° to 45°
- Lower extremity ranges from 0° to 80°
- 110 volt, 50-60Hz motor
- Freight Shipping Specs: 72.5" x 36.6" x 27.6" (184 x 93 x 70cm) and 180lbs (81.5 kgs)



*Hand Controller*



*Caster Control Bar*



*Height Control Bar  
(1 on each side)*

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**





# 3-Section Hi-Lo Table Manual Hydraulic

Function you can use anywhere.

- Hydraulic assist for easy height adjustment.
- Black polyvinyl cover with 2 inch foam.
- Durable, easy to clean and sanitize.
- Removable face plug with contoured opening.
- Durable ivory powder-coated steel frame.
- Heavy-duty casters with individual engagement levers.

3-Section Hi-Lo Table, Manual Hydraulic

NC94700



## Specifications:

- 27" Wide (70cm); 78" (198cm) Length
- Height Range 19" to 36" (48cm to 91cm)
- Weight Limit 500 lbs (225kgs)
- Head section angle ranges from -33° to 85°
- Lower section angle ranges from -24° to 41°
- Freight Shipping specs: 79" x 30" x 21"  
(200 x 75 x 52cm) and 172 lbs (78kg)



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**



# Virtual Showroom

E-COMMERCE AFFILIATE PROGRAM

Quickly & Easily Create an **Online Store** on your Website.



**Virtual Showroom™** is an easy to install, fully-functional e-commerce platform designed specifically for rehab healthcare providers. **Virtual Showroom™** integrates seamlessly into your existing website and gives your clients direct access to the products their therapists have recommended. It is simple to create your own **Virtual Showroom™** and it requires no programming skills.

## *Highlights:*

- Clients stay on your website when shopping for recommended products.
- Customizable: Select Featured Products, Categories, and color-match to your website.
- Adapts to mobile devices and tablets.

## *Benefits:*

- No inventory required.
- Earn 25% commissions on most sales.
- Competitive prices for all products.
- North Coast handles all customer service inquiries.
- Full product warranties and our 100% Satisfaction & Return Policy for all sales.

LEARN MORE



Sign-up at [www.myncmstore.com/virtual-showroom](http://www.myncmstore.com/virtual-showroom)  
or call **800.821.9319**

# RxCart™ 2.0 *Increase Clinic Revenue with a Click.*

With RxCart™ 2.0, therapists populate a shopping cart with the exact products and supplies that are recommended for their client – size, weight, resistance level, orientation and even brand. No more sending patients to the internet jungle, instead, provide a quick and efficient process that will assure your clients get exactly what you want them to have. RxCart™ orders are delivered right to your clients home. No inventory, no hassle, and no burdensome retail transaction.

Harness the power of the therapist recommendation.


Simple to use on a tablet or smartphone while consulting with patients.

Makes it easy for your patient to get the specific therapy products you recommend.

## ***Benefits include:***

- Earn 25% commission on most transactions.
- Track all RxCart™ sales with detailed analytics.
- View status, send reminders for ignored RxCarts.
- View all sent orders, by clinic and therapist.
- No interruptions of your professional practice.
- Full product warranties & NCM's 100% satisfaction & return policy for every sale.

LEARN MORE



Sign-up at [www.mycmstore.com/rx-cart](http://www.mycmstore.com/rx-cart)  
or call **800.821.9319**



# Our Company

For 50 years, North Coast Medical has maintained an outstanding reputation in the global rehabilitation market. We are more than just a wholesale distributor and manufacturer of products, equipment and supplies. North Coast is a strategic supply chain partner, providing innovative opportunities for our customers and business partners.

Now a second-generation family business, we continue to invest in our future while steadfastly adhering to strong ethics and integrity. We invite you to browse our website and discover what is *more than expected* at North Coast Medical.

Thousands of medical and rehabilitative products are available at [www.ncmedical.com](http://www.ncmedical.com) – log in to view your professional pricing. Orders can also be placed by phone at **800.821.9319**, or by fax at **877.213.9300**.

## Proprietary Brands

Throughout the years, we've developed and/or acquired medical devices with trademarks and brand names including:



# North Coast

MEDICAL & REHABILITATION PRODUCTS

*More Than Expected™*