



Your
COMPLETE SYSTEM
for **HAND** and
Upper Extremity
THERAPY

PORTABLE. CONVENIENT.

- Lightweight: ~ 10 lb (4.6 kg) per model
- Easily disinfected
- Ideal for in clinic and Telehealth visits
- Secures to any smooth horizontal surface
- Use seated or standing
- 2 Year Warranty
- Patent Pending

TWO UNIQUE MODELS.
ONE COMPLETE SYSTEM.

STANDARDIZED ACTIVITY SIMULATION

Meet desired outcomes by tracking therapeutic progressions, resistance levels and repetitions.

REAL-LIFE SIMULATION

Address real patient use cases with dynamic 'finger to shoulder' exercises.


DESIGNED BY THERAPISTS

Increase strength, range of motion and functional use with interchangeable, task-specific tools.

IMPROVE FUNCTIONALITY. IMPACT LIVES.

**INTERCHANGEABLE
ATTACHMENTS FACILITATE
OVER 70 SPECIFIC MOVEMENTS.**

To accommodate various therapeutic exercises for targeted muscle groups, FEPSim® is available in both horizontal and vertical orientations. Each includes a unique set of tools to meet designated rehabilitation requirements.



FEPSim 1000 Horizontal
FEP-FS1000H

FEPSim 1000 Vertical
FEP-FS1000V

SPECIFICATIONS:

Standard Attachment Packages:

FEPSim 1000 Horizontal:

- 2-inch Dowel Grips (x2)
- Doorknob
- Lever
- T-handle
- Key
- Storage Base

FEPSim 1000 Vertical:

- Large Lid
- Faucet
- Small Cap
- Screwdriver and Screw
- Wrench and Bolt
- 1.5-inch Dowel Grip
- Storage Base

Indications:

Incorporating therapeutic movements of Flexion, Extension, Pronation and Supination, FEPSim® aids in the rehabilitation of:

- Forearm, Wrist and Hand Fractures
- Neurological Conditions including Stroke
- Tendinopathies of the Upper Extremity
- Upper Extremity Weakness
- Osteoarthritis
- Complex Regional Pain Syndrome
- Pre and Post-Surgical Upper Extremity Conditions

Custom attachments may be ordered.



HORIZONTAL



VERTICAL



Visit www.ncmedical.com
or contact your North Coast Medical
sales representative for details.
Toll-free: 800.821.9319