

THERAPEUTIC  
PAIN RELIEF

# Ortho Thermic™

## Thumb CMC Restriction Splint



## Soothing Comfort

Provides support and warmth for the thumb CMC joint.

Features a moldable metal stay for enhanced support at the thumb's base.

Silicon beads on outer glove surface improve gripping ability.

Heat-generating antimicrobial liner.

### OrthoThermic™ Thumb CMC Restriction Splint

Left	Right	Size	MP Circumference
NC79670	NC79671	Small	6" - 7" (15-18 cm)
NC79672	NC79673	Small +	6½" - 7¼" (17-18 cm)
NC79674	NC79675	Medium	7" - 8" (18-20 cm)
NC79676	NC79677	Medium +	7¼" - 8¼" (18-21 cm)
NC79678	NC79679	Large	8" - 9" (20-23 cm)
NC79680	NC79681	Large +	8¼" - 9½" (21-24 cm)
NC79682	NC79683	X-Large	9" - 10" (23-25 cm)



Visit [www.ncmedical.com](http://www.ncmedical.com) or contact your  
North Coast Medical sales representative for details  
Toll-free: 800.821.9319

COMPRESSION + HEAT