

PATHWAY® CTS 2000

PATHWAY® CTS 1500

PATHWAY® CTS Suite

Pelvic Muscle Rehabilitation Systems



The
PROMETHEUS
G R O U P®

Contact us today for more information.
One Washington Street, Suite 3171, Dover, NH
info@theprogrp.com | theprogrp.com | 1 (800) 442-2325

EXPERIENCE THE DIFFERENCE

The Prometheus Group®, Est. 1989, specializes in the design, development, manufacture, and sales of systems for pelvic health.

At Prometheus, we are passionate about **CUSTOMER VALUE**. With more than 5000 unique clientele worldwide, we believe this is represented in the essence of our business:

Unparalleled Technical Support

We provide the highest level of technical assistance through a number of support channels; phone, email, and remote access. We understand your time is precious. Our extended hours make our support services more accessible; resulting in a positive user experience.

Exceptional Customer Service

Our priority is to build and maintain lasting relationships with our clients. Friendly, knowledgeable, and professional team members are always available to take your call.


Dedicated Sales Managers

Our team of field experts are ready to guide you in the right direction and support you throughout your decision-making process and beyond.

Made in the U.S.A.


The Prometheus Group® combines modern technology and ease of use to provide the most comprehensive and advanced high-quality American-made equipment on the market.

The choice is yours, **EXPERIENCE THE DIFFERENCE**, and become part of the Prometheus family.



PATHWAY® CTS 2000 or PATHWAY® CTS 1500 Pelvic Muscle Rehabilitation Systems Include:

- Notebook Computer
- Color Printer
- USB Serial Interface Cable
- USB Flash Drive Software Program
- PATHWAY® CTS 2000 or CTS 1500 Pelvic Muscle Rehabilitation System Operator's Guide
- PATHWAY® CTS 2000 or CTS 1500 Pelvic Muscle Rehabilitation System Module



PATHWAY® CTS 2000

The Pathway® CTS 2000 Pelvic Muscle Rehabilitation System Module providing one Channel of Vaginal and Rectal Intracavity sEMG or External Pelvic Muscle sEMG, One Channel of Accessory Muscle sEMG, One Channel of Vaginal and Rectal Electrical Stimulation, One Channel of Rectal Pressure Manometry.



PATHWAY® CTS 1500

The Pathway® CTS 1500 Pelvic Muscle Rehabilitation System Module providing one Channel of Vaginal and Rectal Intracavity sEMG or External Pelvic Muscle sEMG, One Channel of Accessory Muscle sEMG, One Channel of Vaginal and Rectal Electrical Stimulation.

Simultaneously Monitor Two Channels of sEMG and Rectal Pressure Manometry

Available on the Pathway® CTS 2000 only.

sEMG and Stimulation on a Single Platform

One channel of sEMG to perform muscle contractions and monitor muscle fatigue while providing synchronous stimulation.

Support Services

Unlimited hardware and software Technical Support and Team Viewer remote support. (On-site clinical training included with the Pathway® CTS 2000 only).

Accessories Included

Starter kits include multiple accessories to begin in clinic use immediately.

Database

Database encryption with automatic or manual backup.

Custom Protocols

Set goals, display type, single or split-screen, colors, work/rest intervals, and markers. Name and store protocols on the treatment list.

Documentation

Chart note and narrative reports. PDF and Microsoft Word editable documentation. Review and compare sessions, track progress, and transfer documentation to a specific file location.

Premier Privacy and Compliance

Multi-level password protection, provider and administrative access controls, unique patient identifiers and audit management.



Software documentation includes an extensive compilation of clinical reports, histories, examinations and exercise prescriptions.



Let's Grow Your Practice Together.

Simultaneous Pelvic and Accessory Muscle sEMG with Rectal Pressure

* Intuitive Software Interface for Seamless Treatment Navigation

Freeze Frame or Restart Study

Comprehensive Menu of User Defined Treatment Parameters

Interactive Event Markers

Drop-Down List of Default and User Defined Protocols

Work/Rest Prompts

Change Ranges on the Fly

Resting, Mean, Minimum and Maximum Pressure Measurements

Real Time Percentage Success

On Screen Audio and Visual Goal Adjustments with Matrixed Audio Feedback

*Rectal pressure available on the PATHWAY® CTS 2000 only.

Automatically Construct Chart Note and Narrative Documents

File View Options Help

Current Patient Visit

Visit #1 on 03/16/2011 at 11:47:57 AM

Attachments:

Intake Medical Issues	Created
History	Not Created
Physical Exams	Not Created
Diagnosis and Plan of Care	Created
Session Evaluation	Not Created
Exercise Prescription	Not Created

Care Administered by: Joe Provider B.S. RN

View/Print/Save Documents

Review Previous Patient Visits

Current Patient Visit

Visit #1 on 03/16/2011 at 11:47:57 AM

Attachments:

Intake Medical Issues	Created
History	Not Created
Physical Exams	Not Created
Diagnosis and Plan of Care	Created
Session Evaluation	Not Created
Exercise Prescription	Not Created

Care Administered by: Joe Provider B.S. RN

View/Print/Save Documents

Review Previous Patient Visits

Track and Optimize Patient Progress

Formulate Home Programs, Add Questionnaires, Histories, and Exams

Synchronize Vaginal or Rectal sEMG with Stimulation

12.5, 50, 100 and 200Hz Treatment Protocols

Immediately Curtail Stimulation

Adjust the Stimulation Amplitude During Treatment

Extensive Library of Treatment Displays to Include Animated Video Clips


Template Emulations

3D Area Scrolls

Bars

Rose Opening and Closing

Space Shuttle Docking into Space Station



Visit Us Online At: THEPROGRP.COM

PATHWAY[®] CTS Suite

Pelvic Muscle Rehabilitation Systems

The Prometheus Group[®]

One Washington Street, Suite 3171

Dover, NH 03820 | Tel: 603-749-0733 | Fax: 603-749-0511

info@theprogrp.com | theprogrp.com



ISO 13485/MDSAP CERTIFIED • ©The Prometheus Group[®]. All rights reserved. MKT-022 Rev. B.
COL, UROGYN, UROL, PT