



ionto⁴

DRUG-DELIVERY IONTOPHORESIS
ELECTRODE TREATMENT KIT

A dependable solution for single-use iontophoresis treatments.

- Made of non-sensitizing, hypoallergenic, medical grade materials.
- Maximum strength, precise buffering with unique carbon-Ag/AgCl (silver/silver chloride) technology.
- Each electrode has a large surface adhesion area with a highly absorbent drug matrix.
- Large ground electrode surface area lowers sensation and increases comfort.
- Box includes 12 individual treatment kits.
- **Kit includes** (1) activity electrode, (1) return electrode and (1) alcohol prep pad.
- Available in 3 Square sizes plus Butterfly.
- Not made with natural rubber latex.



Available in a variety of shapes and sizes:



SMALL
1.5cc
NC89375



MEDIUM
2.5cc
NC89376



LARGE
4.0cc
NC89377



BUTTERFLY
2.0cc
NC89378

Distributed by:

Enter Name

Enter Email

Enter Phone Number