

Norco® 3-Section Hi-Lo Table Manual Hydraulic

Function you can use anywhere.

- Hydraulic assist for easy height adjustment.
- Black polyvinyl matte cover with 2 inch foam.
- Durable, easy to clean and sanitize.
- Removable face plug with contoured opening.
- Durable powder-coated steel construction.
- Heavy-duty casters with individual engagement levers.

NC94700
3-Section Hi Lo Table,
Manual Hydraulic



Hydraulic pump lever

Specifications:

- 27" Wide (70cm); 78" (198cm) Length
- Height Range 19" to 36" (48cm to 91cm)
- Weight Limit 500 lbs (225kgs)
- Head section angle ranges from -33° to 85°
- Lower section angle ranges from -24° to 41°
- Freight Shipping specs: 79" x 30" x 21" (200 x 75 x 52cm) and 172 lbs (78kg)
- Lifetime warranty on steel structures and welds.
3 years warranty on frame, hydraulic pump, upholstery, moveable parts, hardware, and accessories.

Distributed by: