

Norco® 3-Section Hi-Lo Table Electric

Full-feature table for ease of use and client comfort.


- Dual full-length foot bars and hand controller for easy height adjustment.
- Premium matte black polyurethane upholstery offers comfort and durability. Easy to clean and sanitize.
- Superior comfort with 2.5" foam padding. (1.25" standard foam, plus 1.25" Slo-Foam™)
- Removable face plug with contoured opening.
- Convenient paper roll dispenser and cutter for quick setup between clients.
- Durable powder-coated steel construction.
- One-step caster control bar engages all casters at once for easy table repositioning.



NC94710 3-Section Hi-Lo Table, Electric



Paper Dispenser



Caster Control Bar



*Height Control Bar
(1 on each side)*



Hand Controller

Specifications:

- 30" (76cm) Wide; 76" (193cm) Length
- Height Range 19" to 37" (48cm to 94cm)
- Lifting Capacity: 550 lbs (250 kgs)
- Static Weight Capacity: 850 lbs (385 kgs)
- Head section angle ranges from -90° to 45°
- Lower extremity ranges from 0° to 80°
- 110 volt, 50-60Hz motor
- Freight Shipping Specs: 72.5" x 36.6" x 27.6" (184 x 93 x 70cm) and 180lbs (81.5 kgs)
- Lifetime warranty on steel structures and welds.
3 year warranty on motor, upholstery, moveable parts, hardware, and accessories.



Distributed by: **Name**
Email
Phone Number