

Big-Grip™ T-Handle Rocker Knife

Ergonomic grip for easy cutting.

- Adaptive T-Handle Rocker Knife utilizes a rocking motion to help cut food.
- Ideal for those with hand arthritis, neurological impairments, weak grasp, limited hand coordination or generalized arm weakness.
- Features a 1 $\frac{3}{8}$ " (3.5cm) T-shaped handle that increases leverage and control for cutting.
- Soft, built-up handle contains flexible ribbing that adapts to any grip.
- The stainless steel blade measures 4 $\frac{1}{2}$ " (11cm) and can be sharpened with a knife sharpener.
- Includes a protective blade cover to be used while stored.



Non-Slip Soft Grip

NC28207 Big-Grip™ T-Handle Rocker Knife

**DISHWASHER
SAFE
STAINLESS
STEEL**



Not made with
natural rubber latex.

