

Name: \_\_\_\_\_

Phone: \_\_\_\_\_

School/Organization: \_\_\_\_\_

E-mail: \_\_\_\_\_

Dealer/Rep: \_\_\_\_\_

Sales Rep: \_\_\_\_\_ Quote: \_\_\_\_\_

# PT2000 BARIATRIC HI-LO MAT TABLE

### STANDARD FEATURES:

- Lifting capacity: 1000 lbs.
- Height range: 18-33"
- Durable steel base with powder coat finish
- 1.5" High density foam cushion
- Rounded corners
- Single foot control
- Convenient drop down casters
- ADA Compliant



PT2000-47



PT2000-47

### STANDARD TABLE:

**PT2000-47**  
48"W x 84"L x 18-33"H

QTY: \_\_\_\_\_

**PT2000-57**  
60"W x 84"L x 18-33"H

QTY: \_\_\_\_\_

**PT2000-68**  
72"W x 96"L x 18-33"H


QTY: \_\_\_\_\_

### OPTIONS:

**SBO-047** - Side Bumper for Mat Table, 4"W x 7"L

**SBO-057** - Side Bumper for Mat Table, 5"W x 7"L

**SBO-068** - Side Bumper for Mat Table, 6"W x 8"L



### STANDARD UPHOLSTERY CHOICES:

- |  |   |
|--|---|
| <input type="radio"/>  Royal Blue      | <input type="radio"/>  American Beauty |
| <input type="radio"/>  Imperial Blue   | <input type="radio"/>  Lark            |
| <input type="radio"/>  Blue Ridge      | <input type="radio"/>  Burgundy        |
| <input type="radio"/>  Deep Violet     | <input type="radio"/>  Dove            |
| <input type="radio"/>  Forest          | <input type="radio"/>  Charcoal        |
| <input type="radio"/>  Mandarin Orange | <input type="radio"/>  Mocha           |
| <input type="radio"/>  Black           | <input type="radio"/>  Putty           |